

Microscope Slides & Cover Glasses

Conforms to ISO 8037 standard

- Fine ground edges slide are manufactured by state-of-art manufacturing process which results in a high quality and polished ground edges.
- Refractive index of 1.53 ± 0.02
- From extremely uniform, thin and flat glass sheets which are very resistant to corrosion and fogging.
- The slide is used for routine microscopy application.



9100 - MICROSCOPE SLIDES

- Made from the finest white glass
- Uniform quality
- High- Precision ground edges on all four sides
- Thin, Flat and resistant to corrosion or fogging
- Pre-cleaned and ready to use
- Subjected to strict quality control checks, examined individually and environmentally controlled dry before packaging
- Specially designed for meeting the stringent needs of users
- Free from all visual defects

9100P02 - PLAIN



9100D01 - DOUBLE FROSTED



Product Code	Description	Pack Size mm L x B x WT	Quantity Per Packet No.	Packets / Case	Price / Packet ₹
9100P02	Plain, Ground Edges	76 x 26 x 1	50	25	115.00
9100D01	Double Frosted, Ground Edges	76 X 26 X 1	50	25	135.00

L : Length x B : Breadth x WT : Wall Thickness



9115 - COVER GLASS

- Made from the finest white glass with excellent quality
- Uniform quality
- Available in both square and rectangular shape
- Cover glass is cut, cleaned, sorted and weighed
- Resistance to corrosion and fogging
- Free from ragged edges, oily film and finger prints
- Packed in reclosable plastic boxes for protection during storage
- Pre cleaned and ready to use
- Subjected to strict quality control checks, examined individually and environmentally controlled to dry before packaging

9115S01 - SQUARE



9115S02 - SQUARE



9115S03 - SQUARE



Product Code	Description	Pack Size mm	Quantity Per Packet / Piece	Packets / Case	Price / Packet ₹
9115S01	Square	18 x 18	100	50	90.00
9115S02	Square	22 x 22	100	50	90.00
9115S03	Square	24 x 24	100	50	95.00

9115R01 - RECTANGULAR



9115R02 - RECTANGULAR



9115R03 - RECTANGULAR




	Description	Pack Size mm	Quantity Per Packet / Piece	Packets / Case	Price / Packet ₹
9115R01	Rectangular	22 x 40	100	50	95.00
9115R02	Rectangular	22 x 50	100	50	100.00
9115R03	Rectangular	22 x 60	100	50	105.00

BOROSIL®

Lab Accessories

Tongs



Clamp with PVC Coated Prong



Clamp



Clamp with Bosshead



Burette Clamp



Boss Head



Retort Ring & Test Tube Holder



Retort Base



Rod For Retort Base



Z Shape Test Tube Stand



Burner



Spatula

