

BOROSIL[®]

Analytical Vials & Closures

“for Labs where
compliance matters”

Why should I choose BOROSIL vials?

No matter which make of instrument you use, it is critical that your sample is carried in a reliable vial.

Whereas, it is important that vials are made within stringent dimensions in controlled conditions to perform consistently, it is even more critical to a Chemist to be able to confidently rule out the vial as a contributor to an OOS (out of specification) incident.

BOROSIL being a manufacturer of vials, provides a Chemist access to the detailed testing & quality documentation of vials thereby enabling a chemist to investigate root causes effectively.



9mm - SCREW NECK VIALS

USP Type I / Flat Bottom

Product Code	Volume ml	Mouth Dia	Vial Colour	Write on Patch	Via Dia mm	Vial Height	Vial Pack	Vial Case Qty	Price ₹
VO02C00900S000	2	9	Clear	-	12	32	100	10 x 100	10460.00
VO02C10900S000	2	9	Clear	Yes	12	32	100	10 x 100	12510.00
VO02A00900S000	2	9	Amber	-	12	32	100	10 x 100	13820.00
VO02A10900S000	2	9	Amber	Yes	12	32	100	10 x 100	16830.00

Caps for 9mm - 2ml Screw Neck Vials

Screw / PTFE

Product Code	Colour	Cap Dia	Cap MOC	Septa MOC	Septa Amount	Septa Mount / Feature	Vial Pack	Vial Case Qty	Price ₹
CS000009PKS001*	Black	9	PP	Silicone	Friction		100	10 x 100	14010.00
CS000009PBS011	Blue	9	PP	Silicone			100	10 x 100	12140.00
CS000009PBS012	Blue	9	PP	Silicone		Pre-Slit	100	10 x 100	13690.00
CS000009PBS111	Blue	9	PP	Silicone	Bonded		100	10 x 100	14780.00
CS000009PBS112	Blue	9	PP	Silicone		Pre-Slit	100	10 x 100	15560.00
CS000009PBS113	Blue	9	PP	Silicone		Pre-Slit	100	10 x 100	16990.00
CS000009PKS114	Black	9	PP	Polypropylene		Pre-Slit	100	10 x 100	15560.00
CS000009PBS115	Blue	9	PP	Silicone		Pre-Slit	f100	10 x 100	16990.00

* Top - Closed

9mm - 2ml Vial & Cap Kits

Screw / Silicone / PTFE

Product Code	Volume ml	Mouth Dia	Vial Colour	Septa MOC	Septa Mount / Feature	Sets / Pk	Vial Case Qty	Price ₹
VC02C009PKS001*	2	9	Clear	Silicone	Friction	100	10 x 100	26520.00
VC02C009PBS011	2	9	Clear	Silicone		100	10 x 100	22600.00
VC02C009PBS012	2	9	Clear	Silicone	Friction / Pre Slit	100	10 x 100	24150.00
VC02C009PBS111	2	9	Clear	Silicone	Bonded Pre Slit	100	10 x 100	25240.00
VC02C009PBS112	2	9	Clear	Silicone		100	10 x 100	26020.00
VC02C009PBS113	2	9	Clear	Silicone		100	10 x 100	26630.00
VC02C009PKS114	2	9	Clear	Polypropylene		100	10 x 100	26020.00
VC02C009PBS115	2	9	Clear	Silicone		100	10 x 100	27450.00
VC02C109PBS011	2	9	Clear	Silicone	Friction	100	10 x 100	24470.00
VC02C109PKS001*	2	9	Clear	Silicone		100	10 x 100	24650.00
VC02C109PBS012	2	9	Clear	Silicone	Friction / Pre Slit	100	10 x 100	26200.00
VC02C109PBS111	2	9	Clear	Silicone	Bonded	100	10 x 100	27290.00
VC02C109PBS112	2	9	Clear	Silicone	Bonded Pre Slit Write on	100	10 x 100	28070.00
VC02C109PBS113	2	9	Clear	Silicone		100	10 x 100	27450.00
VC02C109PKS114	2	9	Clear	Polypropylene		100	10 x 100	28070.00
VC02C109PBS115	2	9	Clear	Silicone		100	10 x 100	29500.00

* Top - Closed

9mm - 2ml Vial & Cap Kits

Screw / Silicone / PTFE

Product Code	Volume ml	Mouth Dia	Vial Colour	Septa MOC	Septa Mount / Feature	Sets / Pk	Vial Case Qty	Price ₹
VC02A009PKS001*	2	9	Amber	Silicone	Friction	100	10 x 10	27000.00
VC02A009PBS011	2	9	Amber	Silicone		100	10 x 100	25960.00
VC02A009PBS012	2	9	Amber	Silicone	Friction / Pre-Slit	100	10 x 100	27510.00
VC02A009PBS111	2	9	Amber	Silicone	Bonded	100	10 x 100	28600.00
VC02A009PBS112	2	9	Amber	Silicone	Bonded Pre - Slit	100	10 x 100	29380.00
VC02A009PBS113	2	9	Amber	Silicone		100	10 x 100	30810.00
VC02A009PKS114	2	9	Amber	Polypropylene		100	10 x 100	29830.00
VC02A009PBS115	2	9	Amber	Silicone		100	10 x 100	30810.00
VC02A109PKS001*	2	9	Amber	Silicone	Friction	100	10 x 100	27830.00
VC02A109PBS011	2	9	Amber	Silicone		100	10 x 100	28970.00
VC02A109PBS012	2	9	Amber	Silicone	Friction / Pre - Slit	100	10 x 100	30520.00
VC02A109PBS111	2	9	Amber	Silicone	Bonded	100	10 x 100	31610.00
VC02A109PBS112	2	9	Amber	Silicone	Bonded Pre Slit Write on	100	10 x 100	32390.00
VC02A109PBS113	2	9	Amber	Silicone		100	10 x 100	33820.00
VC02A109PKS114	2	9	Amber	Polypropylene		100	10 x 100	32390.00
VC02A109PBS115	2	9	Amber	Silicone		100	10 x 100	33820.00

* Top - Closed



Vials

11mm- 2ml Crimp Neck Vials

USP Type I / Flat Bottom

Product Code	Volume ml	Mouth Dia	Vial Colour	Write on Patch	Via Dia mm	Vial Height	Vial Pack	Vial Case Qty	Price ₹
VO02C01100C000	2	11	Clear	No	12	32	100	10 X 100	9190.00
VO02C11100C000	2	11	Clear	Yes	12	32	100	10 X 100	9930.00
VO02A01100C000	2	11	Amber	No	12	32	100	10 X 100	10770.00
VO02A11100C000	2	11	Amber	Yes	12	32	100	10 X 100	13220.00

Caps for 11mm -2ml Crimp Neck Vials

Silver / Top - Open / Silicone

Product Code	Volume ml	Vial Colour	Cap MOC	Septa MOC	Lining	Septa Mount	Caps Pack	Case Qty	Price ₹
CS000011ASC011	11	Silver	Aluminum	Silicone	PTFE	Friction	100	10 X 100	11200.00
CS000011ASC012*	11	Silver	Aluminum	Silicone	PTFE	Friction	100	10 X 100	12760.00
CS000011ASC021	11	Silver	Aluminum	Silicone	PTFE	Friction	100	10 X 100	8710.00
CS000011ASC022*	11	Silver	Aluminum	Silicone	PTFE	Friction	100	10 X 100	10270.00

* Feature - Pre Slit

11mm Vial and Cap Kits :

Crimp / Open / Silicone / PTFE

Product Code	Volume ml	Neck Dia	Vial Colour	Type	Top	Septa mount	Sets / Pk	Vial Case Qty	Price ₹
VC02C011ASC011	2	11	Clear	Crimp	Open	Friction	100	10 X 100	20390.00
VC02C011ASC012	2	11	Clear	Crimp	Open	Friction	100	10 X 100	21950.00
VC02C011ASC021	2	11	Clear	Crimp	Open	Friction	100	10 X 100	17900.00
VC02C011ASC022	2	11	Clear	Crimp	Open	Friction	100	10 X 100	19460.00
VC02C111ASC011	2	11	Clear	Crimp	Open	Friction	100	10 X 100	21130.00
VC02C111ASC012	2	11	Clear	Crimp	Open	Friction	100	10 X 100	22690.00
VC02C111ASC021	2	11	Clear	Crimp	Open	Friction	100	10 X 100	18640.00
VC02C111ASC022	2	11	Clear	Crimp	Open	Friction	100	10 X 100	20200.00
VC02A011ASC011	2	11	Amber	Crimp	Open	Friction	100	10 X 100	21970.00
VC02A011ASC012	2	11	Amber	Crimp	Open	Friction	100	10 X 100	23530.00
VC02A011ASC021	2	11	Amber	Crimp	Open	Friction	100	10 X 100	19480.00
VC02A011ASC022	2	11	Amber	Crimp	Open	Friction	100	10 X 100	21040.00
VC02A111ASC011	2	11	Amber	Crimp	Open	Friction	100	10 X 100	24420.00
VC02A111ASC012	2	11	Amber	Crimp	Open	Friction	100	10 X 100	25980.00
VC02A111ASC021	2	11	Amber	Crimp	Open	Friction	100	10 X 100	21930.00
VC02A111ASC022	2	11	Amber	Crimp	Open	Friction	100	10 X 100	23490.00

Vials

20mm Headspace GC Vials

USP Type I / Clear /

Product Code	Volume ml	Mouth Dia mm	Neck	Via Dia mm	Vial Height mm	Vial Profile	Vial Pack	Vial Case Qty	Price ₹
VO10F12000C000	10	20	Crimp	22	45	Flat Bottom	100	5 X 100	7220.00
VO10R12000C000	10	20	Crimp	22	45	Round Bottom	100	5 X 100	8600.00
VO20F12000C000	20	20	Crimp	22	75	Flat Bottom	100	5 X 100	14180.00
VO20R12000C000	20	20	Crimp	22	75	Round Bottom	100	5 X 100	14180.00
VO10F12000S000	10	18	Crimp	22	45	Flat Bottom	100	5 X 100	12005.00
VO10R12000S000	10	18	Crimp	22	45	Round Bottom	100	5 X 100	12915.00
VO20F12000S000	20	18	Crimp	22	75	Flat Bottom	100	5 X 100	19180.00
VO20R12000S000	20	18	Crimp	22	75	Round Bottom	100	5 X 100	21165.00

Caps for Headspace GC Vials

Silver Colour / Top - Open / Silicone

Product Code	Cap Dia mm	Cap MOC	Cap Type	Lining	Septa Dia	Feature Height mm	Caps / Pk	Case Qty	Price ₹
CS000020ASC041	20	Aluminium	Crimp	PTFE	Friction	Ultra Clean	100	5 X 100	14810.00
CS000020ASC051	20	Aluminium	Crimp	PTFE	Friction	Everyday	100	5 X 100	16615.00
CS000020MSC041	20	Magnetic	Crimp	PTFE	Friction	Ultra Clean	100	5 X 100	18425.00
CS000020MSS061	18	Magnetic	Screw	PTFE	Friction	Ultra Clean	100	5 X 100	32510.00

Headspace Vial and Cap Kits

Clear / Crimp Type / Top - Open / Silicone / PTFE

Product Code	Volume ml	Neck Dia	Type	Septa Mount	Feature	Sets / Pk	Case Qty	Price ₹
VC20F120ASC041	20	20	Crimp	Friction	Flat Bottom - Ultra Clean	100	5 X 100	30795.00
VC20F120ASC051	20	20	Crimp	Friction	Flat Bottom - Everyday	100	5 X 100	28990.00
VC20F120MSC041	20	20	Crimp	Friction	Flat Bottom - Magnetic	100	5 X 100	32605.00
VC20R120ASC041	20	20	Crimp	Friction	Round Bottom - Ultra Clean	100	5 X 100	30795.00
VC20R120ASC051	20	20	Crimp	Friction	Round Bottom - Everyday	100	5 X 100	28990.00
VC20R120MSC041	20	20	Crimp	Friction	Round Bottom - Magnetic	100	5 X 1003	35980.00
VC10R120ASC041	10	20	Crimp	Friction	Round Bottom - Ultra Clean	100	5 X 100	25215.00
VC10R120ASC051	10	20	Crimp	Friction	Round Bottom - Everyday	100	5 X 100	23410.00
VC10R120MSC041	10	20	Crimp	Friction	Round Bottom - Magnetic	100	5 X 100	27025.00
VC10F120ASC041	10	20	Crimp	Friction	Flat Bottom - Ultra Clean	100	5 X 100	23835.00
VC10F120ASC051	10	20	Crimp	Friction	Flat Bottom - Everyday	100	5 X 100	22030.00
VC10F120MSC041	10	20	Crimp	Friction	Flat Bottom - Magnetic	100	5 X 100	25645.00
VC10F120MSS061	10	18	Screw	Friction	Flat Bottom - Magnetic	100	5 X 100	44515.00
VC20F120MSS061	20	18	Screw	Friction	Flat Bottom - Magnetic	100	5 X 100	51690.00
VC10R120MSS061	10	18	Screw	Friction	Round Bottom - Magnetic	100	5 X 100	45425.00
VC20R120MSS061	20	18	Screw	Friction	Round Bottom - Magnetic	100	5 X 100	53675.00

EPA VIALS



EPA/ TOC Pre-Assembled Vial and Cap Kits

Screw Type / Top - Open / Feature - Pre Assembled

Product Code	Volume ml	Neck Dia mm	Vial Colour	Septa MOC	Lining Mount	Septa Height	Vial Pack mm	Vial Case Qty	ViaPrice ₹
VC20C124PWS111	20	24	Clear	Silicone	PTFE	Bonded	72	1 X 72	8530.00
VC30C124PWS111	30	24	Clear	Silicone	PTFE	Bonded	72	1 X 72	9450.00
VC40C124PWS111	40	24	Clear	Silicone	PTFE	Bonded	72	1 X 72	11380.00
VC60C124PWS111	60	24	Clear	Silicone	PTFE	Bonded	72	1 X 72	15550.00
VC20A124PWS111	20	24	Amber	Silicone	PTFE	Bonded	72	1 X 72	9470.00
VC30A124PWS111	30	24	Amber	Silicone	PTFE	Bonded	72	1 X 72	10590.00
VC40A124PWS111	40	24	Amber	Screw	PTFE	Bonded	72	1 X 72	12960.00
VC60A124PWS111	60	24	Amber	Screw	PTFE	Bonded	72	1 X 72	18040.00
VC40C124POS111*	40	24	Clear	Screw	PTFE	Bonded	72	1 X 72	17150.00

* Pre Assembled & top cover